



TECHNICAL DATASHEET

Item code	7240/INC PROOF	
Name	FROG PLUS - Self loading trolley for incubator predisposed for certified fixing system	
Dimensions [mm]	Length	1950
	Length with folded trolley	1660
	Width	560
	Max height	820
	Folded height	270
	Loading height	600 ÷ 660
Wheels	n°4 Ø200 mm (2 fixed+2 rotating) n°4 Ø125 mm n°2 Ø100 mm	
Weight [kg]	38	
Max load [kg]	250	
Used materials	AISI304, HDPE	
Declared use	Transporting incubators. The trolley is suitable for transporting incubators on vehicles appropriately prepared with fastening systems	
Conformity	MDR 2017/745 – Class I Medical Device EN 1789 certified	



2306

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it