



TECHNICAL DATASHEET

| | |
|------------------------|---|
| Item code | 7240 PROOF |
| Name | FROG PLUS – Self-loading stretcher with Trendelenburg certified |
| Dimensions [mm] | Length 1950 |
| | Length with folded trolley 1660 |
| | Width 560 |
| | Max height 920 |
| | Folded height 325 |
| | Loading height 600 ÷ 660 |
| Weight [kg] | 44 |
| Max load [kg] | 250 |
| Used materials | AISI304, ABS, HDPE |
| Wheels | n°4 Ø200 mm (2 fixed+2 rotating) n°4 Ø125 mm n°2 Ø100 mm |
| Declared use | Transport, in the lying down position, of temporarily or permanently disabled patients with pathologies that allow this. The stretcher is suitable for transporting patients on vehicles appropriately prepared with fastening system |
| Conformity | MDR 2017/745 – Class I Medical Device EN 1865-1 – EN 1789 certified |



2306

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it