



TECHNICAL DATASHEET

| | | |
|------------------------|--|-----------|
| Item code | 7260 PROOF | |
| Name | FROG LITE - Self loading aluminium stretcher with trendelenburg certified | |
| Dimensions [mm] | Length | 1960 |
| | Length with folded trolley | 1670 |
| | Width | 560 |
| | Max height | 915 |
| | Folded height | 310 |
| | Loading height | 600 ÷ 660 |
| Weight [kg] | 35,1 | |
| Max load [kg] | 250 | |
| Used materials | Al alloys, HDPE | |
| Wheels | n°4 Ø200 mm (2 fixed+2 rotating) n°4 Ø125 mm n°2 Ø100 mm | |
| Declared use | Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled | |
| Conformity | MDR 2017/745 – Class I Medical Device EN 1865-1 – EN 1789 certified | |



2306

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it