



TECHNICAL DATASHEET

Item code	7265 PROOF	
Name	FROG LITE ADVANCE - Aluminium self loading stretcher with extractable handles and trendelenburg	
Dimensions [mm]	Length	1960
	Length with folded trolley	1670
	Width	560
	Max height	915
	Folded height	310
	Loading height	600 ÷ 660
Weight [kg]	36,5	
Max load [kg]	250	
Used materials	Al alloys, HDPE	
Wheels	n°4 Ø200 mm (2 fixed w/brake +2 rotating) n°4 Ø125 mm n°2 Ø100 mm	
Declared use	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled	
Conformity	MDR 2017/745 – Class I Medical Device EN 1865-1 – EN 1789 certified	



2306

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it