



TECHNICAL DATASHEET

Item code	7270 PROOF
Name	MBS-01 MEBER BARIATRIC STRETCHER - Self loading stretcher for bariatric transport
Dimensions [mm]	Length 1970
	Width 560
	Max height 920
	Folded height 250
	Loading height 600
Weight [kg]	47,5
Max load [kg]	350
Used materials	AISI 304, ABS, HDPE
Wheels	n°4 Ø200 mm (2 fixed w/brake +2 rotating) n°2 Ø125 mm n°2 Ø100 mm
Declared use	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled
Conformity	MDR 2017/745 – Class I Medical Device EN 1865-3 – EN 1789 certified



2306

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it