



TECHNICAL DATASHEET

Item code	7300	
Name	MEBER-X – Multilevel stretcher	
Dimensions [mm]	Length	1660 ÷ 1970
	Width	600
	Height	405 · 555 · 670 · 775 · 870 · 965 · 1060 · 1145 · 1230
	Loading height	0 · 170 · 290 · 385 · 480 · 580 · 670 · 760 · 840
Weight [kg]	45,9	
Max load [kg]	320	
Used materials	Al alloys, polyamide	
Characteristics	▪ 4 rotating wheels ▪ Rear wheels with brake ▪ Moving wheels: nr. 4 Ø200 ▪ Loading wheels: nr 2 Ø100 ▪ Adjustable backrest: 0 ÷ 80° ▪ Trendelenburg: 15° ▪ Fowler: 35°	
Declared use	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled	
Conformity	MDR 2017/745 – Class I Medical Device	



2306

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it