



TECHNICAL DATASHEET

Item code	910/GIN PROOF	
Name	WINNER – Self-loading trolley for incubators predisposed for certified fixing system	
Dimensions [mm]	Length	1940
	Length with folded trolley	1550
	Width	560
	Max height	800
	Folded height	220
	Loading height	620 ÷ 650
Used materials	Al alloys, S235JR, ABS, HDPE	
Wheels	n°4 Ø150 mm (2 fixed w/brake + 2 rotating) n°2 Ø125 mm n°2 Ø100 mm	
Weight [kg]	31	
Max load [kg]	170	
Declared use	Transport of incubators and loading on rescue vehicles	
Conformity	MDR 2017/745 – Class I Medical Device EN 1789 certified	



2207

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it