



## TECHNICAL DATASHEET

<b>Item code</b>	<b>910 PROOF</b>	
<b>Name</b>	<b>WINNER – Certified self-loading stretcher</b>	
<b>Dimensions [mm]</b>	Length	1940
	Length with folded trolley	1550
	Width	560
	Max height	875
	Folded height	260
	Loading height	620 ÷ 650
<b>Used materials</b>	Al alloys, S235JR, ABS, HDPE	
<b>Wheels</b>	n°4 Ø150 mm (2 fixed w/brake + 2 rotating) n°2 Ø125 mm n°2 Ø100 mm	
<b>Weight [kg]</b>	37	
<b>Max load [kg]</b>	170	
<b>Declared use</b>	Transport in a lying position and loading on rescue vehicle of patients with appropriate conditions, who are temporarily or permanently disabled	
<b>Conformity</b>	MDR 2017/745 – Class I Medical Device EN 1865-1 – EN 1789 certified	



2207

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:  
[info@meber.it](mailto:info@meber.it)