



TECHNICAL DATASHEET

Item code	910
Name	WINNER – Self-loading stretcher
Dimensions [mm]	Length 1940
	Length with folded trolley 1620
	Width 560
	Max height 875
	Folded height 260
	Loading height 450 ÷ 650
Used materials	Al alloys, S235JR, ABS, HDPE
Wheels	n°4 Ø150 mm (2 fixed w/brake + 2 rotating) n°2 Ø125 mm
Weight [kg]	34
Max load [kg]	170
Declared use	Transporting of incubators and loading on rescue vehicles
Conformity	MDR 2017/745 – Class I Medical Device



2207

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it