



## TECHNICAL DATASHEET

|                        |  |           |
|------------------------|--|-----------|
| <b>Item code</b>       | <b>914 PROOF</b>   |           |
| <b>Name</b>            | <b>WINNER – Certified self-loading stretcher with Trendelenburg and fowler</b>   |           |
| <b>Dimensions [mm]</b> | Length   | 1940      |
|                        | Length with folded trolley   | 1550      |
|                        | Width  | 560       |
|                        | Max height   | 875       |
|                        | Folded height  | 260       |
|                        | Loading height   | 620 ÷ 650 |
| <b>Used materials</b>  | Al alloys, S235JR, ABS, HDPE   |           |
| <b>Wheels</b>          | n°4 Ø150 mm (2 fixed w/brake + 2 rotating)<br>n°2 Ø125 mm<br>n°2 Ø100 mm   |           |
| <b>Weight [kg]</b>     | 40   |           |
| <b>Max load [kg]</b>   | 170  |           |
| <b>Declared use</b>    | Transport in a lying position and loading on rescue vehicle of patients with appropriate conditions, who are temporarily or permanently disabled |           |
| <b>Conformity</b>      | MDR 2017/745 – Class I Medical Device<br>EN 1865-1 – EN 1789 certified   |           |



2207

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:  
[info@meber.it](mailto:info@meber.it)