



TECHNICAL DATASHEET

Item code	914	
Name	WINNER – Self-loading stretcher with Trendelenburg and fowler	
Dimensions [mm]	Length	1940
	Length with folded trolley	1620
	Width	560
	Max height	875
	Folded height	260
	Loading height	450 ÷ 650
Used materials	Al alloys, S235JR, ABS, HDPE	
Wheels	n°4 Ø150 mm (2 fixed w/brake + 2 rotating) n°2 Ø125 mm n°2 Ø100 mm	
Weight [kg]	37	
Max load [kg]	170	
Declared use	Transport in a lying position and loading on rescue vehicle of patients with appropriate conditions, who are temporarily or permanently disabled	
Conformity	MDR 2017/745 – Class I Medical Device	



2207

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it