

KUMA

POWERED STAIR CHAIR

Art. 668

Integrated light to illuminate the stairs
in low light conditions

Command
console



Art. 668

KUMA

Black powered stair chair

Kuma is the entry level model in MeBer's line of powered stair chairs. Frame in aluminium, lithium-ion battery powered, front handgrips, ergonomic armrests, telescopic grab bar (with a small console for track speed control) and 4 castors (2 fixed rear Ø200 mm with brakes and 2 front swivels castors Ø125 mm). The electric motor operates the tracks, allowing an effortless ascent and descent of the stairs. Kuma has a 181 kg load capacity.

Supplied with battery, battery charger, seat cushion, headrest, chest belt, ventral belt. Practical and easy to manoeuvre, the KUMA is extremely easy to use. Colour: black. **Device compliant with MDR 2017/745 - Class I DM.**

181 kg
MAXIMUM
LOADING

30+
AUTONOMY
OVER
30 RAMPS

KUMA

POWERED STAIR CHAIR



Foldable handles

Telescopic front handles



Headrest

Command console

Foldable handles

Motorized tracks



All rights are reserved. Me.Ber. reserves the right to modify the content without prior notice. Me.Ber. srl Unipersonale cannot be held responsible for any damage caused by incorrect or incomplete information. Pictures shown are for illustration purpose only.

Accessories

Art. 12180

Battery for KUMA Powered Stair Chair



Art. 12182

Battery charger for KUMA Powered Stair Chair



Art. 977

Fastening system for KUMA Powered Stair Chair



Supplied with

Battery
Charger
Seat cushion
Headrest
Belts

Art.	Material	Length	Width	Height	Seat height	Weight	Loading capacity	Battery	Folded thickness	Autonomy
668	Aluminium alloy	79 cm	52 cm	108 - 164 cm	50 cm	34 kg	181 kg	24 V	24 cm	over 30 flights of 20 steps with a load of 130 kg



Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270 43124 Fontanini (Parma) Italy
Ph. +39 0521 648770 | Fax +39 0521 648780 | e-mail: info@meber.it

www.meber.it