



## TECHNICAL DATASHEET

|                        |  |           |
|------------------------|--|-----------|
| <b>Item code</b>       | <b>7070/4RG PROOF</b>  |           |
| <b>Name</b>            | <b>MERCURY - Stainless steel certified variable height self-loading stretcher with 4 rotating wheels</b>             |           |
| <b>Dimensions [mm]</b> | Length   | 1970      |
|                        | Length with folded trolley   | 1650      |
|                        | Width  | 570       |
|                        | Max height   | 950       |
|                        | Intermediate height  | 720       |
|                        | Folded height  | 310       |
|                        | Loading height   | 640 ÷ 660 |
| <b>Used materials</b>  | AISI 304, ABS, HDPE  |           |
| <b>Wheels</b>          | n°4 Ø200 mm<br>n°2 Ø150 mm<br>n°2 Ø60 mm   |           |
| <b>Weight [kg]</b>     | 50,6   |           |
| <b>Max load [kg]</b>   | 250  |           |
| <b>Declared use</b>    | Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled |           |
| <b>Conformity</b>      | MDR 2017/745 - Class I Medical Device<br>EN 1865-1 - EN 1789 certified   |           |



2604

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes.

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:  
[info@meber.it](mailto:info@meber.it)