



TECHNICAL DATASHEET

Item code	679	
Name	EXTRA LITE ZENA - Reduced width chair, 4 foldable wheels with footrest, armrest and aluminium platform	
Dimensions [mm]	Length	650 ÷ 1220
	Width	470 (rear levers closed) 500 (rear levers open)
	Height	1020
	Folded thickness	290
	Seat height	560
Weight [kg]	14,2	
Max load [kg]	230	
Used materials	Al alloys	
Wheels	n°2 Ø200 mm (fixed) n°2 Ø125 mm (rotating w/brake)	
Declared use	Manual transportation of a seated patient with use-compatible pathologies who are temporarily or permanently disabled. The product is NOT fit for use inside rescue vehicles with patient on board	
Conformity	MDR 2017/745 – Class I Medical Device EN 1789 Certified	



2605

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it