



TECHNICAL DATASHEET

Item code	7077/4RG PROOF
Name	MERCURY 5 - Stainless steel certified high version 5 position self loading stretcher with 4 rotating wheels
Dimensions [mm]	Length 1970
	Length with folded trolley 1670
	Width 570
	Height 1010 – 910 – 820 – 700 – 390
	Folded height 310
	Loading height 690 ÷ 710
Used materials	AISI 304, aluminium, HDPE
Wheels	n°4 Ø200 mm n°2 Ø150 mm n°2 Ø60 mm
Weight [kg]	49,8
Max load [kg]	250
Declared use	Transporting, in the lying down position, of temporarily or permanently disabled patients with pathologies that allow this. The stretcher is suitable for transporting patients on vehicles appropriately prepared with fastening systems
Conformity	MDR2017/745 - Class I Medical Device EN 1865-1 - EN 1789 certified



2605

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes.

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it