



## TECHNICAL DATASHEET

<b>Item code</b>	<b>7092 PROOF</b>																
<b>Name</b>	<b>MERCURY EL - Black painted steel cot certified EN 1865</b>																
<b>Dimensions [mm]</b>	<table border="1"> <tr> <td>A</td> <td>1970</td> <td>A1</td> <td>-</td> </tr> <tr> <td>B</td> <td>570</td> <td></td> <td></td> </tr> <tr> <td>C</td> <td colspan="3">950 - 310 (folded)</td> </tr> <tr> <td>D</td> <td colspan="3">640 ÷ 660</td> </tr> </table>	A	1970	A1	-	B	570			C	950 - 310 (folded)			D	640 ÷ 660		
A	1970	A1	-														
B	570																
C	950 - 310 (folded)																
D	640 ÷ 660																
<b>Weight [kg]</b>	46,5																
<b>Max load [kg]</b>	250																
<b>Used materials</b>	S235JR, ABS, HDPE																
<b>Wheels</b>	n°4 Ø200 mm (2 fixed+2 swivel) n°2 Ø150 mm n°2 Ø60 mm																
<b>Features</b>	Trendelenburg: 15° tilt																
<b>Declared use</b>	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled																
<b>Conformity</b>	MDR 2017/745 – Class I Medical Device EN 1865-1 - EN 1789 certified																



2510

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

**Me.Ber. s.r.l. Unipersonale**

Strada Langhirano, 270  
43124 Fontanini (Parma) Italy  
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:  
[info@meber.it](mailto:info@meber.it)