



TECHNICAL DATASHEET

Item code	7094 PROOF																
Name	MERCURY EL - Black painted steel cot high version certified EN 1865																
Dimensions [mm]	<table border="1"> <tr> <td>A</td> <td>1970</td> <td>A1</td> <td>-</td> </tr> <tr> <td>B</td> <td>570</td> <td></td> <td></td> </tr> <tr> <td>C</td> <td colspan="3">990 · 310 (folded)</td> </tr> <tr> <td>D</td> <td colspan="3">690 ÷ 710</td> </tr> </table>	A	1970	A1	-	B	570			C	990 · 310 (folded)			D	690 ÷ 710		
A	1970	A1	-														
B	570																
C	990 · 310 (folded)																
D	690 ÷ 710																
Weight [kg]	46,6																
Max load [kg]	250																
Used materials	S235JR, ABS, HDPE																
Wheels	n°4 Ø200 mm (2 fixed+2 swivel) n°2 Ø150 mm n°2 Ø60 mm																
Features	Trendelenburg: 15° tilt																
Declared use	Transport, in a prone position, of patients with appropriate conditions, who are temporarily or permanently disabled																
Conformity	MDR 2017/745 – Class I Medical Device EN 1865-1 - EN 1789 certified																



2510

All rights are reserved. Any variations on images and subjects can be done without notice. Me.Ber. s.r.l. Unipersonale is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes

Me.Ber. s.r.l. Unipersonale

Strada Langhirano, 270
43124 Fontanini (Parma) Italy
Tel. +39 0521 648770 | Fax +39 0521 648780 | e-mail:
info@meber.it