



7260 proof Frog Lite 7260 PROOF

"FROG LITE" self loading aluminium stretcher with trendelenburg certified



FEATURES

Made of painted aluminium alloy, extremely light, it is also provided with trendelenburg. The 4 Ø 200 mm polyurethane wheels (2 fixed with brake and 2 rotating with self-centring system) grant high comfort on uneven surfaces. Easy to clean and desinfect.

- · Rigid platform made of painted aluminium alloy
- Reclining backrest 0°-80°, operated by a gas piston, lockable in any position
- Folding sides
- · Vertical legs protected by plastic buffets
- · Rotating wheels are self centring
- · Collapsable frontal trolley

In coupling with the proper fixation Art. 955 Proof, it is 10g resistant. Supplied complete with mattress (Art. 7013) and belts (nr. 2 Art. 608-A/MEB o 608-N/MEB and nr. 1 Art. 612/AN-2/MEB or Art. 612/N-2/MEB).



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

MAIN FEATURES

10G Certification	₩
Height (cm)	91,5
Max loading height (cm)	66
Min loading height (cm)	60
Folded height (cm)	31
Width (cm)	56
Length (cm)	196
Weight (kg)	36,9
Loading capacity (kg)	250
Filter feature	
Material	Painted aluminium HDPE (High-density polyethylene)

IMAGES





SEE ALSO



7765 Proof Aluminium self loading stretcher "Frog lite Advance" with exctractable handles and Trendelenburg



Frog 7240 PROOF
"FROG PLUS" self loading stretcher
with trendelenburg certified



Frog 7210 PROOF "FROG" yellow self loading stretcher with trendelenburg certified

ACCESSORIES



955 PROOF EN 1789 fixing system for Proof rollin stretchers



Instrument support for stretchers



936
Phleboclysis stand for automatic loading stretcher Winner and Frog



12025Black multi usage object holder for automatic loading stretcher



12106Mattress for AIR-COM stretchers Anti-bedsore mattress with
alternate segments



7015Folding sky mattress black for bariatric MBP platform



12154"BEAR PODS" Bariatric containment system for self-loading stretchers



707Red 6-handle fireproof cover with
Me.Ber tape for self-loading
stretcher



938 Oximeb oxygen cylinder holder



12054 Oximeb Maxi - Oxygen cylinder holder



12128Adjustable system for headrest of Meber-X and Mercury stretchers



12080
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



12080-001Monting kit for bariatric platform on self loading stretcher



12070/BG "COVER OX 2-3 BC"- professional cylinder holder for 2 and 3 lt oxygen bottles



12070/BR COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



12072/BG
"COVER OX 5 BG"- professional
cylinder holder for 5 lt oxygen
bottles



12072/BR
"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



12084
Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel

STANDARD ACCESSORIES



7013/N 3 pcs. Skay fire-proof electrowelded mattress



612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



608/N-MEB 2 pieces black belt with metalic hook and MeBer band (L. 155 cm)