



11010/LN Cormos

Portable aspirator for secretions with disposable liner



FEATURES

Medical aspirator for secretions CORMOS, compact and powerful, indicated for emergency services. CORMOS works directly with the 12V. cc. of the vehicle or with an inner accumulator with more than 45 min. of autonomy. CORMOS is silent and powerful thanks to its engine of 25 litres/min which enables a constant aspiration, with control manometer, positioned on the front of the aspirator, superior to 850 mbar. Easy to carry thanks to its particular and compact dimension, the anatomic handle and its limited weight of 4,6 Kg. New case design, made of ABS/UL high resistance anti-shock with dimension of cm 34x11x24. It is supplied with 1000 ml disposable liner, green MeBer bag, suction rigid probe with hose and rechargeable bracket 10C.

The "CORMOS" aspirator for secretions is manufactured according to the ISO 10079-1:2009 product normative, IIb class, EN 1789 and is certified CE 0123, according to the directive 93/42/CEE for medical devices and european directive 104 (automotive).

MAIN FEATURES

Width (cm)	34
Height (cm)	24
Depth (cm.)	11
Weight (kg)	4,6

SEE ALSO



Cormos Plus 11022/FA White complete medical suction unit with autoclavable canister



Cormos Plus 11022/LN White complete medical suction unit with canister for disposable unit



Cormos 11010/FA
Portable aspirator for secretions
with autoclavable canister

ACCESSORIES



11162 Disposable 1-litre canister



11024 Serres disposable liners 1000 ml (x 6 pz/pcs.)



11025 Serres disposable liners 1000 ml (x 12 pz./pcs.)



11026 Serres disposable liners 1000 ml (x 36 pz/pcs.)



11012 Autoclavable canister



11013 Antibacterical filter



11034 Battery charger for Cormos 2 poles



11022-003 2 poles cable for aspirator CORMOS

STANDARD ACCESSORIES



11014 Suction rigid probe and hose for "CORMOS" and "CORMOS PLUS" aspirators



11016Green carrying bag for suction unit



11028 Autoclavable canister



11023
10g certified wall bracket with integrated charging system for CORMOS and CORMOS PLUS