



## 9234 Artis

Cervical collar set 4 pieces with carrying case.  
Art. 9220 - 9224 - 9226 - 9228



### FEATURES

Rigid cervical collar in 1 piece made with stamped polyethylene coloured to identify various sizes and E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening and central hole for tracheotomy. Available in 6 different sizes.  
Art. 9220 - 9224 - 9226 - 9228

### MAIN FEATURES

Trauma



### SEE ALSO



**Artis 9220**  
Rigid cervical collar in 1 piece size  
INFANT



**Artis 9224**  
Rigid cervical collar in 1 piece size  
NECKLESS



**Artis 9226**  
Rigid cervical collar in 1 piece size  
SHORT



**Artis 9228**  
Rigid cervical collar in 1 piece size  
REGULAR