



910/GIN Winner

Self loading trolley for incubators



FEATURES

Designed to load and transport incubators on rescue vehicles. This trolley is part of "Winner" range, keeping the all same features. Made of anodized aluminium.

MAIN FEATURES

Pediatric	₩
Height (cm)	87,5
Max loading height (cm)	65
Min loading height (cm)	45
Folded height (cm)	26
Width (cm)	56
Length (cm)	194
Length with folded trolley (cm)	162
Weight (kg)	28
Material	Anodized aluminium
Loading capacity (kg)	170

SEE ALSO



Winner 910/GIN PROOF Self loading trolley for incubators predisposed for certified fixing system

ACCESSORIES



Rear fixing system Rear fixing system for roll-in H=148 mm



Front blocking device Front blocking device for automatic loading stretcher



12048 Certified interface for incubator fixation for self loading trolleys



12050 Extended trolley incubator interface