




## 9040 Galaxy-P


Yellow spinal board with pins



### FEATURES

Same features as “Galaxy” spinal boards Art. 9010, 9012 and 9015 equipped with 10 plastic PINS useful for a fast fixing of restraint systems with karabiners. Supplied with 3 belts (605/MEB).

### MAIN FEATURES

Length (cm)	183
Width (cm)	41
Thickness (cm)	5
Loading capacity (kg)	150
Weight (kg)	5,8
Trauma	
Usage	Immobilization
Material	ABS

---

ACCESSORIES

---



**Jelly 629**

Orange standard head immobilizer  
for spinal board



**574**

Blue bag with 3 belts art. 572

---

STANDARD ACCESSORIES

---



**605/MEB**

50 mm orange belt in one piece  
with plastic buckle and ME.BER.  
ribbon