



873 Sirius

Cervical collar set 3 pieces with carrying case
Art. 864 - 866 - 868



FEATURES

Rigid cervical collar in 1 piece made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening which different colours indicate the various sizes. Central hole for tracheotomy. Available in 6 different sizes.
Art. 864 - 866 - 868

MAIN FEATURES

| | |
|--------|----------------|
| Usage | Immobilization |
| Trauma | |

SEE ALSO



Sirius 876
Cervical collar set 6 pieces with carrying case



Antares 856
Cervical collar set 6 pieces with carrying case



Antares 854
Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



Sirius 864
Rigid cervical collar in 1 piece size STOUT



Sirius 866
Rigid cervical collar in 1 piece size STANDARD



Sirius 868
Rigid cervical collar in 1 piece size MEDIUM
