



19116

180 mm ramp for stretcher supports



FEATURES

180 mm ramp with raised sides.

MAIN FEATURES

Length (cm)	18
Width (cm)	63
Height (cm)	1,7
Weight (kg)	2,3

ACCESSORIES



19130
"MTX" single space stretcher
support with lateral movement
with ramp mm. 180 - EN 1789
certified



19130-955
"MTX-955" single space support
with lateral movement and ramp of
mm. 180, equipped with EN1789
certified fixing system art. 955 EN1789 certified



19130-956

"MTX-956" single space support
with lateral movement and ramp of
mm. 180, equipped with EN1789
certified fixing system art. 956 EN1789 certified



19136 "MX" fixed single space stretcher support with ramp of mm 180 -EN1789 certified



19136-955
"MX-955" fixed single space stretcher support with ramp of mm. 180, equipped with certified fixing system art. 955 - EN1789 certified



19136-956
"MX-956" fixed single space stretcher support with ramp of mm. 180, equipped with en1789 certified fixing system art. 956 - EN1789 certified