



19112

100 mm ramp for stretcher supports



FEATURES

100 mm ramp with raised sides.

MAIN FEATURES

Length (cm)	10
Width (cm)	62
Height (cm)	5
Weight (kg)	1,4

ACCESSORY OF



19130

"MTX" single space stretcher support with lateral movement with ramp mm. 180 - EN 1789 certified



19130-955

"MTX-955" single space support with lateral movement and ramp of mm. 180, equipped with EN1789 certified fixing system art. 955 - EN1789 certified



19130-956

"MTX-956" single space support with lateral movement and ramp of mm. 180, equipped with EN1789 certified fixing system art. 956 - EN1789 certified



19136

"MX" fixed single space stretcher support with ramp of mm 180 - EN1789 certified



19136-955

"MX-955" fixed single space stretcher support with ramp of mm. 180, equipped with certified fixing system art. 955 - EN1789 certified



19136-956

"MX-956" fixed single space stretcher support with ramp of mm. 180, equipped with en1789 certified fixing system art. 956 - EN1789 certified