



7015

Folding sky mattress black for bariatric MBP platform



FEATURES

Mattress for bariatric MBP platform made of expanded polyurethane and fireproof, waterproof and anti-tear skay cover. It is foldable into 2 parts for an easy storage and it is supplied with 6 Velcro strips for anchoring it to the platform frame.

Dimensions: Lenght x Width x Thickness: 173 x 96 x 5,5 (max thickness 10) cm

*All Pictures shown are for illustration purpose only.

MAIN FEATURES

Bariatric	₩
Length (cm)	173
Width (cm)	96
Thickness (cm)	5,5 10

IMAGES



SEE ALSO



12080-001Monting kit for bariatric platform on self loading stretcher

ACCESSORY OF



12080
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



MBS-01 MEBER BARIATRIC STRETCHER 7270 Self loading stretcher for bariatric transport