



975

“Ergolift fixing system” - Fixing system for chair Extra Ergolift



FEATURES

Fixing device designed to ensure the fixation of Extra Ergolift chair inside the ambulance patient compartment. It is made of stainless steel and it is corrosion resistant and easy to clean. The innovative monobloc configuration allows the operator to lock and unlock the chair fastly and easily.

The fixing system for Extra Ergolift chair can be also used as fixing system for the traditional Extra chair: 670/BR, 676 and 678.

Device certified according to the European harmonized standards UNI EN 1789.



MAIN FEATURES

Width (cm)	58
Height (cm)	111
Depth (cm.)	23
Weight (kg)	8

IMAGES



SEE ALSO



960
"Extra fixing system" - fixing system
for chair Extra

ACCESSORY OF



Extra Ergolift 671
Powered stair chair



Extra 670/BR
Track assisted stair chair with arms



Extra Lite Standard 676
Foldable carrying chair with 4
wheels, footrest and arms



678 Extra Lite Basic
Extra Lite Basic, foldable carrying
chair
