



## 891/SG


Flake basic independent chambers vacuum mattress



### FEATURES

Flake basic independent chambers vacuum mattress 200x90.

### MAIN FEATURES

Length (cm)	200
Width (cm)	90
Weight (kg)	6,2
Loading capacity (kg)	150
Trauma	
Size	Adults
Material	PVC-coated polyester

### SEE ALSO



**"Flake" 891**  
Independent chambers vacuum mattress



**Chesterfield 9114**  
Vacuum mattress with independent chambers with central welding