



# 8396-001

Pedestal for art. 8396 – 8398



## FEATURES

Main characteristics: yellow cabinet pedestalfor art. 8396-8398 in carbon steel P02 7/10 thick powder coated in epoxy polyester with holes for fixing to the floor.

## MAIN FEATURES

Height (cm)	80
Depth (cm.)	120
Width (cm)	27

## SEE ALSO



8432-001  
Pedestal for art. 8432

---

ACCESSORY OF

---



**8398**

Defibrillator metal cabinet with thermoregulation and alarm for outer use



**8396**

Defibrillator metal cabinet with thermoregulation and alarm for outer use

---