



8540

Multiparametric monitor TG



FEATURES

The TG monitor is a professional tool with many operational options and memory, fitted with accessories and circuits.

The monitor has a dual function: diagnosis and monitoring.

Allows for the individual or simultaneous detection of various vital parameters.


5 parameters

- 5 channel, 2 track ECG
- Systolic/diastolic/average blood pressure
- Heartbeat
- Pulse oximeter
- Breathing control
- Twin probe temperature

Functions

- Multi-parameter monitor with number and track display.
- Operates on battery and mains supply - Power supply: 100 ÷ 240 V ac, 50 ÷ 60 Hz. It is also equipped with a battery to allow continuity of use without external power supply (up to 4 hours with the standard battery)
- Complete with built-in printer
- Network connection for remote management
- Parameter adjustment and setting
- 10" monitor

MAIN FEATURES

| Cardio |  |
|-------------|---|
| Height (cm) | 24,6 |
| Width (cm) | 26,1 |
| Depth (cm.) | 14,6 |
| Weight (kg) | 2,8 |

ACCESSORIES



8540-001
Reusable Spo2 sensor for adults for
multi-parameter monitor



8540-002
27-35 cm adult cuff



8540-005
Rechargeable lithium battery



8540-006
Thermal graph paper



8540-003
NIBP bracelet, E7, pediatric, 16-21.5
cm



8540-004
NIBP cuff, E5, neonatal, 10-15 cm

STANDARD ACCESSORIES



8540-001
Reusable Spo2 sensor for adults for
multi-parameter monitor



8540-002
27-35 cm adult cuff