



# 1319

Emergency bag “Giulietta”



## FEATURES

Red Emergency bag made in closed-cell PVC, with PVC bottom, plastic supports and upper transparent pocket for customizations.

Compact but roomy, thanks to its internal modularity it allows the organization of space in a personalized way in order to optimally house what is necessary for emergencies. 4 pockets, 2 on the sides, one front and one rear, predisposed with supports and rubber bands.

It is also equipped with an additional removable bag with transparent front. Thanks to its construction, it allows internal cleaning and disinfection. Portable by hand, shoulder strap or as a backpack.

## MAIN FEATURES

|             |     |
|-------------|-----|
| Width (cm)  | 45  |
| Height (cm) | 33  |
| Depth (cm.) | 34  |
| Weight (kg) | 2,2 |

---

IMAGES

---



---

SEE ALSO

---



**1318**  
Rescue backpack "Romeo"

---