



935

Instrument support for X-frame and Meber Max stretchers



FEATURES

Support for medical equipment such as defibrillator or monitor to be anchored to the X-frame cot's chassis. Supplied with straps for the instalment of medical equipment. Foldable for a better storage.

MAIN FEATURES

Height (cm)	40
Depth (cm.)	34
Epaisseur pliée (cm)	14
Max Width (cm)	63
Min Width (cm)	56

ACCESSORY OF



MeBer-X 7300
"MeBer-X" Multi-level stretcher
complete with mattress and belts