



429

Splint for pediatric limbs



FEATURES

Rigid splints in fire-proof and water-resistant PVC. This device ensures a correct limb immobilization. Equipped with 4 VELCRO straps (different lengths), ensuring a safe and quick application. This model is shaped for a firm lacing to ankles and upper leg and is equipped with integrated foot support. Radiolucent.

Available sizes:

Adult arm / pediatric leg: length 52 cm – base 9 cm Adult leg: length 74 cm – base 13 cm

Pediatric limb: lenght 30 cm, base 9 cm

MAIN FEATURES

Pediatric	₩
Trauma	₩
Material	PVC
Length (cm)	30
Usage	Immobilization

SEE ALSO



Arcos long rigid splint for lower limbs



427 Arcos short rigid splint for upper limbs



431 Arcos rigid splint set with carrying bag