



8506

Pulse oximeter Oxy-110



FEATURES

- Large 3.5 inch auto rotation TFT display for SpO2, PI, Plethysmogram, PR and temperature patented MoveOxy™ SpO2 technology with artifact removal and anti-motion performance
- various choice of SpO2 probes for adult, pediatric and neonate
- built-in lithium battery up to 18 hours continuous monitoring
- adjustable alarm limit, pitch tone level, key tone and brightness
- up to 580 hours SpO2 data storage
- upgradeable to infrared ear temperature

Supplied with:

- cover
- adult clip probe
- charging base
- USB cable

MAIN FEATURES

Thickness (cm)	2,5
Height (cm)	15,8
Width (cm)	7,3
Weight (kg)	0,23
Cardio	₩

IMAGES



SEE ALSO



11220 Oxytest

OXYTEST portable pulse oximeter complete with adult sensor, USB cable and AAA battery



11200 Vital Test Handheld pulse oxymeter complete with adult sensor