



8548

Multiparameter monitor PC-300 SpO2, NIPB



FEATURES

PC-300 handheld Spot-check monitor measures SpO2, Pulse Rate, NIBP and TEMP. Widely used in hospitals and clinics. Equipped with Bluetooth and USB interface, allows to transfer data to PC, tablets and smartphone (Android and iOS).

- up to 6000 groups NIBP data storage, up to 100 different patients (each user can store 999 measurement records)
- quick ear probe TEMP measurement result in 5 seconds
- data transfer to PC/Smartphone by USB or Bluetooth
- optional NIBP pediatric cuff and SpO2 pediatric probe
- optional ECG and Glucometer
- built-in rechargeable lithium battery
- communication protocols available on request

Power Supply

Battery: 3.7 V built-in Li-ion battery Operating voltage: 3.6-4.2 VDC Operating current: =800 mA

Charger: 100 - 240 VAC, 50/60 Hz, 5 VDC or USB

Supplied with:

SpO2 adult probe NIBP adult cuff Ear thermometer (probe covers not included) Lithium battery with charger Carrying bag User manual

MAIN FEATURES

Cardio	₩
Length (cm)	16,5
Thickness (cm)	6,8
Width (cm)	9.6

IMAGES



SEE ALSO



B1 CLASSIC 8534 Automatic blood pressure with adult cuff