




19140-955

"MTX-2-955" double space stretcher support with lateral movement and ramp of mm. 350, equipped with fixing system art. 955 - EN1789 certified



FEATURES

Double space stretcher support with lateral movement and ramp mm 350 in the underneath space can be stored both a scoop stretcher and a spinal board. On request supplied with UNI EN 1789 certified blocking devices Art. 955 or 956. Gas spring included.



Lateral movement: cm 30

MAIN FEATURES

Side sliding	
Width (cm)	63
Length (cm)	229
Weight (kg)	84,5
Loading capacity (kg)	400
Height (cm)	27,5

IMAGES



ACCESSORIES



19110
Stainless steel plate mm 1940 x 572
for stretcher supports



787L
Gas spring for ramp



19123
Ribaltina 130 mm con attacchi per
supporti e rullo



19115
Ribaltina da 180 mm con attacchi
per supporti e rullo

STANDARD ACCESSORIES



821G
Gas spring for ramp