



618

Fixing system for sedan chairs



FEATURES

Designed to fix automatically the sedan chairs mod.650-653-658-659 inside the ambulance. Composed by 2 pieces. Made with stainless steel.



Device certified according to the European harmonized standards UNI EN 1789.

Also available separately 618/INF - 618/SUP

MAIN FEATURES

Length (cm)	46
Depth (cm.)	4
Height (cm)	2,5
Weight (kg)	2,3

SEE ALSO



Maya 650
Folding sedan chair



Maya 653
Folding sedan chair with 4 wheels,
footboards and arms

ACCESSORIES



Quick 659/Y
Folding sedan chair with 4 wheels,
footboards and arms - yellow/black



Quick 658/Y
Folding sedan chair with 2 wheels –
yellow/black
