



## 624 Double

Head immobilizer for scoop stretcher



### FEATURES

This head immobilizer can be mounted and separated on the scoop stretcher. It is supplied with the standard side rigid blocks. Supplied with padded belts for the fixing. It is produced with expanded material at closed cells with an inside aluminium sheet and it is covered by a liquid-proof and bacterial-proof material.

### MAIN FEATURES

Size	Adults
Usage	Immobilization
Width (cm)	40
Depth (cm.)	25
Height (cm)	15
Weight (kg)	0,7
Trauma	

### IMAGES



---

ACCESSORY OF

---



**Maxima 631**  
Grey anodized anti trauma scoop  
stretcher with head/trunk in 1 pc



**Maxima 630**  
Grey anodized anti trauma scoop  
stretcher

---