



630 Maxima

Grey anodized anti trauma scoop stretcher



FEATURES

Thanks to its features the scoop stretcher can be placed under the patient with ease, without lifting the body. Built with carefully designed materials, MAXIMA ensures solidity, reliability and durability.

- Efficient splitting system with locking device
- Adjustable length with automatic locking device
- Bag with n. 3 straps (quick release) for the immobilization of the patient (603/MEB)
- Easily arranged inside the ambulance
- Built with light and sturdy aluminium
- Carrying case available on request (Art. 630/CS)



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

MAIN FEATURES

Size	Adults
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	167
Max length (cm)	202
Folded length (cm)	120
Epaisseur pliée (cm)	9
Weight (kg)	10
Loading capacity (kg)	170
10G Certification	₩
Material	Anodized aluminium
Trauma	₩
Alpine	₩

IMAGES





SEE ALSO



Maxima 630/N Black painted anti trauma scoop



Maxima 630/BBlue painted anti trauma scoop stretcher



Maxima 630/AOrange painted anti trauma scoop stretcher



Maxima 631 Grey anodized anti trauma scoop stretcher with head/trunk in 1 pc

ACCESSORIES









Double 624 Head immobilizer for scoop stretcher



630/CSCarrying case for scoop stretchers



603/MEB

Nylon bag containing 3 belts Art.
607/Ameb

ACCESSORY OF



603/MEBNylon bag containing 3 belts Art.
607/Ameb