



630/A Maxima

Orange painted anti trauma scoop stretcher



FEATURES

Keeping the very same features of Art.630, thanks to a special paint (designed for optimal durability), this stretcher can be used to differentiate operators or a specific area of operation.



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

- Efficient splitting system with locking device
- Adjustable length with automatic locking device
- \bullet Bag with n. 3 straps (quick release) for the immobilization of the patient (603/MEB)
- Easily arranged inside the ambulance
- Built with light and sturdy aluminium
- Carrying case available on request (Art. 630/CS)

MAIN FEATURES

Size	Adults
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	167
Max length (cm)	202
Folded length (cm)	120
Epaisseur pliée (cm)	9
Weight (kg)	10
Loading capacity (kg)	170
10G Certification	✓
Material	Painted aluminium

SEE ALSO



Maxima 631 Grey anodized anti trauma scoop stretcher with head/trunk in 1 pc



Maxima 630/NBlack painted anti trauma scoop stretcher



Maxima 630 Grey anodized anti trauma scoop stretcher



Maxima 630/B Blue painted anti trauma scoop stretcher

ACCESSORIES



616Certified vertical fastening system for scoop stretchers



616/ACertified vertical fixing system for scoop stretchers



614Certified horizontal fastening system for scoop stretchers



630/CSCarrying case for scoop stretchers



Double 624 Head immobilizer for scoop stretcher



603/MEBNylon bag containing 3 belts Art.
607/Ameb