



# 630/B Maxima

---

Blue painted anti trauma scoop stretcher



---

## FEATURES

Keeping the very same features of Art.630, thanks to a special paint (designed for optimal durability), this stretcher can be used to differentiate operators or a specific area of operation.




Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

- Efficient splitting system with locking device
- Adjustable length with automatic locking device
- Bag with n. 3 straps (quick release) for the immobilization of the patient (603/MEB)
- Easily arranged inside the ambulance
- Built with light and sturdy aluminium
- Carrying case available on request (Art. 630/CS)

---

## MAIN FEATURES

---

|                       |   |
|-----------------------|---|
| Size                  | Adults  |
| Usage                 | Rescue  |
| Height (cm)           | 7   |
| Width (cm)            | 42  |
| Length (cm)           | 167   |
| Max length (cm)       | 202   |
| Epaisseur pliée (cm)  | 9   |
| Weight (kg)           | 10  |
| Loading capacity (kg) | 170   |
| 10G Certification     |  |
| Material              | Painted aluminium   |

---

## SEE ALSO

---



**Maxima 631**  
Grey anodized anti trauma scoop  
stretcher with head/trunk in 1 pc



**Maxima 630/N**  
Black painted anti trauma scoop  
stretcher



**Maxima 630**  
Grey anodized anti trauma scoop  
stretcher



**Maxima 630/A**  
Orange painted anti trauma scoop  
stretcher

---

## ACCESSORIES

---



**616/A**  
Certified vertical fixing system for  
scoop stretchers



**630/CS**  
Carrying case for scoop stretchers



**614**  
Certified horizontal fastening  
system for scoop stretchers



**616**  
Certified vertical fastening system  
for scoop stretchers



**Double 624**  
Head immobilizer for scoop  
stretcher



**603/MEB**  
Nylon bag containing 3 belts Art.  
607/Ameb