



## 646 Folding stretcher


Collapsible rescue stretcher foldable in 2 pieces



### FEATURES

Designed for the transport also of heavy patients, this stretcher can be folded in two pieces. The canvas is built in fire-resistant, impermeable, tear-resistant PVC-coated nylon. Fixed handles. Straps are not included. (see "Restraint systems" Art. 604/MEB).

### MAIN FEATURES

Size	Adults
Usage	Transport
Foldable	
Min height (cm)	14
Width (cm)	50
Max length (cm)	201
Folded length (cm)	101
Epaisseur pliée (cm)	17
Weight (kg)	6,2
Loading capacity (kg)	170
Material	Anodized aluminium Nylon

### IMAGES



---

## SEE ALSO

---



**Folding stretcher 643**  
Collapsible rescue stretcher  
foldable in 4 pieces

---

---

## ACCESSORIES

---



**604/MEB**  
Orange bag with 2 straps with  
plastic buckle and ME.BER. ribbon

---