



653 Maya

Folding sedan chair with 4 wheels, footboards and arms



FEATURES

Designed to help rescue personnel during rescue procedures. Made of anodized aluminium. Supplied with 4 wheels (2 rotating with brake Ø 100 mm and 2 fixed Ø 125 mm), easy and safe open/close system, 4 rear handles and 2 front extractable ones, back rest and seat made of thermoplastic material easy to clean and footrest. Supplied complete with black bag containing 3 belts (Art. 590/N).



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

MAIN FEATURES

Foldable	
Footboards and arms	
Wheel type	2 Rotating Wheels + 2 Small Fixed Wheels
Type	
Height (cm)	90
Width (cm)	52
Length (cm)	75
Max length (cm)	131
Epaisseur pliée (cm)	24
Weight (kg)	14
Loading capacity (kg)	170
Folded height (cm)	98
10G Certification	
Cardio	
Material	Anodized aluminium

IMAGES



SEE ALSO



Maya 650
Folding sedan chair

ACCESSORIES



618
Fixing system for sedan chairs

STANDARD ACCESSORIES



590/N

Black bag with 3 belts for sedan
chairs
