



## 653 Maya

Folding sedan chair with 4 wheels, footboards and arms



### FEATURES

Designed to help rescue personnel during rescue procedures. Made of anodized aluminium. Supplied with 4 wheels (2 rotating with brake Ø 100 mm and 2 fixed Ø 125 mm), easy and safe open/close system, 4 rear handles and 2 front extractable ones, back rest and seat made of thermoplastic material easy to clean and footrest. Supplied complete with black bag containing 3 belts (Art. 590/N).



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

### MAIN FEATURES

10G Certification	
Foldable	
Footboards and arms	
Wheel type	2 Rotating Wheels + 2 Small Fixed Wheels
Type	
Height (cm)	90
Width (cm)	52
Length (cm)	75
Max length (cm)	131
Folded thickness (cm)	24
Weight (kg)	14
Loading capacity (kg)	170
Folded height (cm)	98
Material	Anodized aluminium
Seat length (cm)	43
Cardio	

---

## IMAGES

---



---

## SEE ALSO

---



**Maya 650**  
Folding sedan chair

---

## ACCESSORIES

---



**618**  
Fixing system for sedan chairs

---

---

## STANDARD ACCESSORIES

---



**590/N**

Black bag with 3 belts for sedan  
chairs

---