



7072 PROOF Mercury

Certified stainless steel variable height stretcher high version



FEATURES

Designed to take different configurations, it is provided with a multi level system using 2 sliding pistons on each axle. Fixing of stretcher inside the rescue vehicle is granted by certified system Art. 955 Proof (see "Fixing Systems"). Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

Features:


- Ergonomic controls (for rescue vehicle loading operations and for configuration adjustments of the patient support surface)
- Adjustable backrest with lockable gas spring for inclinations from 0° to 75°
- Fully folding 600 mm sides
- HDPE buffers for loading onto the rescue vehicle (on both front and back legs)
- 2 fixed front wheels, Ø 200 mm with parking brakes
- 2 rear castors, 360° Ø 200 mm self centring
- 4 loading wheels of which 2 are Ø 100 mm and 2 are Ø 125 mm (enhanced loading system)
- Equipped with mattress Art. 7013 and belts (2 Art. 608-A/MEB o 608-N/MEB and 1 Art. 612/AN-2/MEB or Art. 612/N-2/MEB)

This version has the same features of the standard model with the only difference that it is possible to load the stretcher on ambulances with higher loading height.



Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

MAIN FEATURES

Size	Adults
Usage	Transport
10G Certification	
Multilevel	
Lying surface	With Trendelenburg
Wheel type	2 Rotating Wheels + 2 Big Fixed Wheels
Type	Self loading stretcher
Intermediate height (cm)	75
Max height (cm)	99
Max loading height (cm)	71
Min loading height (cm)	69
Folded height (cm)	31
Width (cm)	57
Length (cm)	197
Length with folded trolley (cm)	165
Weight (kg)	43
Loading capacity (kg)	250
Material	Stainless steel ABS HDPE (High-density polyethylene)

IMAGES



SEE ALSO



Mercury 7070/4RG PROOF
Stainless steel variable height
stretcher certified with 4 rotating
wheels



Mercury EL 7092 PROOF
Certified self-loading stretcher



Mercury EL 7094 PROOF
Certified self-loading stretcher -
high version



Mercury 7072/4RG PROOF
Certified stainless steel variable
height stretcher high version with
4 rotating wheels



Mercury 7070 PROOF
Certified stainless steel variable
height stretcher

ACCESSORIES



955 PROOF
EN 1789 fixing system for Proof roll-in stretchers



7014/A
4 pcs. Skay fire-proof electrowelded mattress



939
Instrument support for stretchers



938
Oximeb oxygen cylinder holder



12080
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



12080-001
Monting kit for bariatric platform on self loading stretcher



12018N
Telescopic foldable phlebotoclysis stand for self loading stretchers



12154
"BEAR PODS" Bariatric containment system for self-loading stretchers



12070/BG
"COVER OX 2-3 BG" - professional cylinder holder for 2 and 3 lt oxygen bottles



12070/BR
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



12072/BG
"COVER OX 5 BG" - professional cylinder holder for 5 lt oxygen bottles



12072/BR
"COVER OX 5 BR" - professional cylinder holder for 5 lt oxygen bottles



12106
Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



7013/A
3 pcs. Skay fire-proof electrowelded mattress



7013/N
3 pcs. Skay fire-proof electrowelded mattress



12054
Oximeb Maxi - Oxygen cylinder holder



707
Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



12084-001
Cover for bariatric wings set



12170
Black mattress for right bariatric wing for self-loading stretchers



12172
Black mattress for left bariatric wing for self-loading stretchers



608/A-MEB

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



608/N-MEB

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



612/AN-2/MEB

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



12025

Black multi usage object holder for automatic loading stretcher



12200

Bariatric wings set for Mercury self-loading stretchers complete with mattresses, belts with reel and assembly kit

STANDARD ACCESSORIES



12025

Black multi usage object holder for automatic loading stretcher



7013/A

3 pcs. Skay fire-proof electrowelded mattress



612/AN-2/MEB

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



608/A-MEB

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)