



7210 PROOF Frog

“FROG” yellow self loading stretcher with trendelenburg certified



FEATURES

Made of welded stainless steel powder coated, this is a really strong and reliable device. The Ø 200 mm wheels (2 fixed with brake and 2 rotating with self-centring system) are made of soft material that allows high comfort level to the patient. This model of stretcher has trendelenburg position to 0° from 20° too.

- Rigid platform made of reinforced thermoplastic material. · Reclining backrest 0°- 80°, operated by a gas piston, can be locked in any position.
- Folding sides.
- Vertical legs protected by plastic buffers.
- Casters are self-centring so that the stretcher can be positioned in the ambulance without them having to be turned manually.
- Front trolley can be lowered to reduce overall dimensions of the stretcher.

In coupling with the proper fixation Art. 955 Proof, it is 10g resistant.

This stretcher is supplied complete with mattress (Art. 7013) and belts (nr. 2 Art. 608-A/MEB o 608-N/MEB and nr. 1 Art. 612/AN-2 o Art. 612/N-2).



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789.

MAIN FEATURES

| | |
|---------------------------------|---|
| Size | Adults |
| Usage | Transport |
| 10G Certification |  |
| Lying surface | With Trendelenburg |
| Wheel type | 2 Rotating Wheels + 2 Big Fixed Wheels |
| Type | Self loading stretcher |
| Height (cm) | 92 |
| Max loading height (cm) | 66 |
| Min loading height (cm) | 60 |
| Folded height (cm) | 32,5 |
| Width (cm) | 56 |
| Length (cm) | 195 |
| Length with folded trolley (cm) | 166 |
| Weight (kg) | 47,5 |
| Loading capacity (kg) | 250 |
| Material | Painted steel |
| Filter feature | |

SEE ALSO



7265 Proof
Aluminium self loading stretcher
"Frog lite Advance" with
extractable handles and
Trendelenburg



Frog Lite 7260 PROOF
"FROG LITE" self loading aluminium
stretcher with trendelenburg
certified



Frog 7210
"FROG" yellow self loading stretcher
with trendelenburg



Frog 7240 PROOF
"FROG PLUS" self loading stretcher
with trendelenburg certified

ACCESSORIES



955 PROOF
EN 1789 fixing system for Proof roll-in stretchers



936
Phleboclysis stand for automatic loading stretcher Winner and Frog



939
Instrument support for stretchers



938
Oximeb oxygen cylinder holder



12025
Black multi usage object holder for automatic loading stretcher



612/AN-2/MEB
Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



612/N-2/MEB
Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



608/A-MEB
2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



608/N-MEB
2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



12084
Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel



7013/N
3 pcs. Skay fire-proof electrowelded mattress



12154
"BEAR PODS" Bariatric containment system for self-loading stretchers



12106
Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



12018N
Telescopic foldable phleboclysis stand for self loading stretchers



12080
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



12080-001
Monting kit for bariatric platform on self loading stretcher



12070/BG
"COVER OX 2-3 BG" - professional cylinder holder for 2 and 3 lt oxygen bottles



12070/BR
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



12072/BG
"COVER OX 5 BG"- professional
cylinder holder for 5 lt oxygen
bottles



12072/BR
"COVER OX 5 BR"- professional
cylinder holder for 5 lt oxygen
bottles

STANDARD ACCESSORIES



7013/N
3 pcs. Skay fire-proof electrowelded
mattress



612/AN-2/MEB
Adjustable thoracic abdominal belt
orange/black with MeBer band and
4 points fastening for self loading
stretchers and sedan chairs



608/A-MEB
2 pieces orange belt with metallic
hook and MeBer band (L. 155 cm)