



## 7240 PROOF Frog

---

“FROG PLUS” self loading stretcher with trendelenburg certified



---

### FEATURES

Made with electro polished STAINLESS STEEL. Supplied with 4 wheels Ø 200 mm (2 fixed with brake and 2 rotating self centring). The inclination of the legs is designed to give the stretcher better stability during use. The trendelenburg position is operated by a gas piston and can be locked in any position from 0° to 15°.

- Rigid platform made of reinforced thermoplastic material.
- Reclining backrest 0°- 80°, operated by a gas piston, can be locked in any position.
- Folding sides.
- Vertical legs protected by plastic buffers.
- Casters are self-centring so that the stretcher can be positioned in the ambulance without them having to be turned manually.
- Front trolley can be lowered to reduce overall dimensions of the stretcher.

In coupling with the proper fixation Art. 955 Proof, it is 10G resistant.

Supplied complete with mattress (Art. 7013) and belts (nr. 2 Art. 608-A/MEB o 608-N/MEB and nr. 1 Art. 612/AN-2/MEB or Art. 612/N-2/MEB)




Device certified according to the European harmonized standards UNI EN1865-UNI EN1789

---

## MAIN FEATURES

---

Usage	Transport
10G Certification	
Lying surface	With Trendelenburg
Wheel type	2 Rotating Wheels + 2 Big Fixed Wheels
Type	
Height (cm)	92
Max loading height (cm)	66
Folded height (cm)	32,5
Min loading height (cm)	60
Width (cm)	56
Length (cm)	195
Length with folded trolley (cm)	166
Weight (kg)	47,5
Loading capacity (kg)	250
Material	Stainless steel
Filter feature	

---

## SEE ALSO

---



**Frog Lite 7260 PROOF**  
"FROG LITE" self loading aluminium  
stretcher with trendelenburg  
certified



**Frog 7210 PROOF**  
"FROG" yellow self loading stretcher  
with trendelenburg certified



**Frog 7240/INC PROOF**  
"FROG PLUS INC PROOF" self  
loading trolley for incubators  
predisposed for certified fixing  
system



**7265 Proof**  
Aluminium self loading stretcher  
"Frog lite Advance" with  
extractable handles and  
Trendelenburg

## ACCESSORIES



**936**

Phlebotomy stand for automatic loading stretcher Winner and Frog



**939**

Instrument support for stretchers



**938**

Oximeb oxygen cylinder holder



**12080**

"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



**955 PROOF**

EN 1789 fixing system for Proof roll-in stretchers



**12080-001**

Monting kit for bariatric platform on self loading stretcher



**12025**

Black multi usage object holder for automatic loading stretcher



**612/AN-2/MEB**

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



**612/N-2/MEB**

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



**608/A-MEB**

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



**608/N-MEB**

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



**7013/N**

3 pcs. Skay fire-proof electrowelded mattress



**12084**

Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel



**12154**

"BEAR PODS" Bariatric containment system for self-loading stretchers



**12106**

Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



**12018N**

Telescopic foldable phlebotomy stand for self loading stretchers



**12070/BG**

"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



**12070/BR**

COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**12072/BG**  
"COVER OX 5 BG"- professional  
cylinder holder for 5 lt oxygen  
bottles



**12072/BR**  
"COVER OX 5 BR"- professional  
cylinder holder for 5 lt oxygen  
bottles

---

## STANDARD ACCESSORIES

---



**7013/N**  
3 pcs. Skay fire-proof electrowelded  
mattress



**612/AN-2/MEB**  
Adjustable thoracic abdominal belt  
orange/black with MeBer band and  
4 points fastening for self loading  
stretchers and sedan chairs



**608/A-MEB**  
2 pieces orange belt with metallic  
hook and MeBer band (L. 155 cm)