



8014

Compact aneroid sphigmomanometer



FEATURES

The handmade sphigmomanometer has high quality finish. Its components are made applying the highest technological solutions: German movement follows the highest quality standards. Polycarbonate dial glass, alcohol and mechanical strain proof. Sphigmomanometer drain valve easy to activate and very sensitive. Natural rubber bladder. Cuff made by a strong anti-allergic cotton. Manometer, fixed on the cuff, is highly resistant to falls or casual blows. Made in Italy according to the EN 1060-1 and EN 1060-2 rules.

MAIN FEATURES

Size

Adults

SEE ALSO



8008
One hand aneroid
sphigmomanometer



8004
Colored aneroid
sphigmomanometer with carrying
case