



8020/G

Yellow flat phonendoscope



FEATURES

Flat phonendoscope with anodized stamped and lathed aluminium head diameter 46 mm. Threaded metal ring for membrane blocking. Metal chromium plated binaural, length 150 mm. Y rubber tube in one piece, diameter 5 mm length 60 cm. Single printed carton box. Available in 3 different colors.

MAIN FEATURES

| | |
|--------|--------|
| Size | Adults |
| Cardio | |

SEE ALSO



8030/B
Blue stethophonendoscope
"Rappaport"



8025/N
Black stethophonendoscope



8020/N
Black flat phonendoscope



8025/G
Yellow stethophonendoscope



8020/R
Red flat phonendoscope



8025/R
Red stethophonendoscope