



844 Antares


Rigid cervical collar in 2 pieces size STOUT



FEATURES

ANTARES cervical collar is an highly professional instrument indispensable for the emergency operator. It is composed by 2 pieces to be put on the front and back part of the neck enabling thus the perfect immobilization of the rachis. Central hole for tracheotomy and easy positioning thanks to the rip-off fastening which different colours indicate the various sizes. It is made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different sizes.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	
Max length (cm)	36,5
Length (cm)	36
Max height (cm)	19
Min height (cm)	14


SEE ALSO



Antares 856
Cervical collar set 6 pieces with carrying case



Sirius 876
Cervical collar set 6 pieces with carrying case



Sirius 874
Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868



Antares 854
Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



Antares 853
Cervical collar set 3 pieces with carrying case Art. 844 - 846 - 848